

Section Seven

Elastics & Power Products


- page 2** Short Sticks, Molded “O”s & Chopped “O”s
- page 3** Power Chain, Tieback Modules, Power Thread,
Power Tube & Separators
- page 4** Zoo Pack Elastics
- page 6** Zoo Pack Dispenser

Short Sticks

Short Sticks are the newest members of the Molded "O" family. Short Sticks accommodate single-patient use to avoid cross contamination between patients.

- Size .120" available in 24 vibrant colors
- Coordinates with Colored Power Chain
- Flat cross section for improved tie-wing fit
- Package of 1,000 (100 trees/10 per tree)

Purple	640-1240
Orange	640-1241
Blue	640-1242
Pink	640-1243
Black	640-1244
Red	640-1245
Neon Green	640-1246
Yellow	640-1247
White	640-1248
Tooth	640-1249
Teal	640-1250
Rose	640-1251
Violet	640-1252
Emerald Green	640-1255
Maroon	640-1256
Gold	640-1257
Silver	640-1258
Jade Green	640-1259
Ruby Red	640-1260
Cobalt Blue	640-1261
Navy	640-1262
Smoke	640-1263
Clear	640-1264
Gray	640-1265
Regular Sampler Pack	640-1270
Metallica Sampler Pack	640-1271



Short Sticks Sectional Storage System and Shade Sheet

Store your Ormco Short Sticks conveniently with the Short Sticks Sectional Storage System. It contains 24 sections, enough to store each color of Ormco Short Sticks, and comes complete with a double-sided, perforated Shade Sheet with which to label each section.

Sectional Storage System with Shade Sheet
640-1275

Additional Shade Sheet
070-0529



Long regarded as the industry standard of orthodontic elastomers, Ormco's latex free elastomers incorporate a life span dating signature that ensures the reliability of use. Manufacturing our own elastomeric products allows for greater quality control, longer shelf life, and better clinical results for you.

Molded "O"s


- Flat cross section for improved tie-wing fit
- Package of 1,000 (20 trees/50 per tree)

Size	Clear	Gray
.110"	640-0074	640-0077
.120"	640-0075	640-0078


Colored Molded "O"s

- Size .120" available in 17 vibrant colors
- Coordinates with Colored Power Chain
- Flat cross section for improved tie-wing fit
- Package of 1,000 (20 trees/50 per tree)

Purple	640-0100
Orange	640-0101
Blue	640-0102
Pink	640-0103
Black	640-0104
Red	640-0105
Neon Green	640-0106
Yellow	640-0107
White	640-0108
Teal	640-0112
Violet	640-0114
Emerald Green	640-0118
Navy	640-0125
Metallica Silver	640-0121
Metallica Jade Green	640-0122
Metallica Cobalt Blue	640-0124
Light Blue	640-0155



Chopped "O"s



- Most economical elastic ligature
- Easy-dispensing preloaded canes
- Package of 1,000 (20 canes of 50)

Size	Clear	Gray
.090"	637-1092	637-1090
.110"	637-2112	637-2110
.120"	637-3122	637-3120

Power Chain

Original Power Chain

Power Chain I and II are easy-to-use continuous chains of modules for ligation, consolidation and anterior retraction.

- Roll of 15 ft
- Open and wide spacing

	Clear	Gray
Power Chain I (wide space)	639-0012	639-0010
Power Chain II (open space)	639-0013	639-0011

Colored Power Chain

Colored Power Chain features a unique urethane formulation that provides constant application of light forces.

- 19 vibrant colors
- Coordinate with Short Sticks colors
- Open and closed spacing
- Roll of 13.1 ft



Color	Open	Closed
Pink	639-0030	639-0040
Purple	639-0031	639-0041
Black	639-0032	639-0042
Green	639-0033	639-0043
Red	639-0034	639-0044
Blue	639-0035	639-0045
Orange	639-0036	639-0046
Yellow	639-0037	639-0047
White	639-0038	639-0048
Teal	639-0050	639-0060
Rose	639-0051	639-0061
Violet	639-0052	639-0062
Emerald Green	639-0053	639-0063
Silver	639-0054	639-0064
Maroon	639-0055	639-0065
Gold	639-0056	639-0066
Jade Green	639-0057	639-0067
Ruby Red	639-0058	639-0068
Cobalt Blue	639-0059	639-0069

Generation II Power Chain

This advanced elastomeric provides a more uniform continuous force over a longer period of time.

- Closed, open, and wide spacing
- Clinical leader in chain elastics
- Roll of 13.1 ft

	Clear	Gray
Generation II (closed space)	639-0002	639-0001
Generation II (open space)	639-0004	639-0003
Generation II (wide space)	639-0019	639-0018

Tieback Modules

Damon Tieback Modules

Designed exclusively for the Damon™ System, the Damon Tieback Modules are manufactured with the same fatigue-resistant material as Ormco's popular Generation II Chain.

- Optimized force-memory and elongation properties provide consistent uniform forces throughout treatment
- Tear-resistant
- Available in two sizes
- Available in two colors (gray and clear) to complement Damon 2, Damon 3 and Damon 3MX bracket systems
- Package of 40

Length	Clear	Gray
6.1 mm	637-0012	637-0011
8.6 mm	637-0022	637-0021



Power Thread & Tube

Power Thread™

- Consistent, reliable performance
- Ideal for rotation, corrections and cuspid retractions

Size	Length	Clear	Gray
.025"	25 ft	640-1026	640-1025
.030"	25 ft	640-1031	640-1030
.025"	50 ft	640-2026	640-2025
.030"	50 ft	640-2031	640-2030



Power Tube™


Hollow channel allows for secure knotting.

Size	Length	Clear	Gray
.025"	25 ft	640-4026	640-4025
.030"	25 ft	640-4031	640-4030
.025"	50 ft	640-6026	640-6025
.030"	50 ft	640-6031	640-6030



Separators

- Radiopaque
- Provide gentle, continuous forces
- Effective in both anterior and posterior sections
- Highly visible blue color



Pack of 1,000 640-0080

Zoo Pack Elastics

- High-quality surgical latex
- Consistent force pull
- Force ratings based on elastic being extended three times listed diameter

- Also available in 10 sizes of assorted Color Paks (pink, orange, green and purple)
- Zoo packs contain 50 bags of 100





































Shaded sizes are available in Ormco Color Paks

Chipmunk	636-0113	Bear	636-0122
Rabbit	636-0114	Impala	636-0125
Fox	636-0115	Ram	636-0126
Penguin	636-0116	Moose	636-0127
Monkey	636-0117	Buffalo	636-0128

Color Paks contain 40 packs of 100


Intraoral Elastics

Extraoral Elastics

Size	Light	Medium	Med./Heavy	Heavy	Ex./Heavy	Light	Heavy
Strength	2 oz/60 g	3 oz/85 g	3.5 oz/100 g	4.5 oz/130 g	6 oz/170 g	8 oz/230 g	14 oz/400 g
1/8" 3.18 mm	 Hummingbird 630-0010		 Chipmunk 630-0030				
3/16" 4.76 mm	 Quail 630-0011	 Otter 630-0020	 Rabbit 630-0031	 Kangaroo 630-0040	 Impala 630-0050	 Cougar 635-0058	
1/4" 6.35 mm	 Owl 630-0012	 Seal 630-0021	 Fox 630-0032	 Bear 630-0041	 Ram 630-0051	 Leopard 635-0059	
5/16" 7.94 mm	 Parrot 630-0013	 Dolphin 630-0022	 Penguin 630-0033	 Zebra 630-0042	 Moose 630-0052	 Panther 635-0060	 Walrus 635-0065
3/8" 9.35 mm	 Road Runner 630-0014	 Turtle 630-0023	 Monkey 630-0034	 Camel 630-0043	 Buffalo 630-0053	 Tiger 635-0061	 Elephant 635-0066
1/2" 12.7 mm	 Peacock 630-0015		 Donkey 630-0035			 Lion 635-0062	 Whale 635-0067
5/8" 15.9 mm	 Eagle 630-0016		 Llama 630-0036				
3/4" 19.1 mm	 Ostrich 630-0017		 Giraffe 630-0037				



Animal Zoo Packs	Pouches of 15 Zoo Packs	Master Pouches of 50 Zoo Packs	Color Packs of 40 Zoo Packs
Light – 2 oz			
Hummingbird		630-0010	
Quail	635-0031	630-0011	
Owl	635-0032	630-0012	
Parrot		630-0013	
Road Runner		630-0014	
Peacock		630-0015	
Eagle		630-0016	
Ostrich	635-0037	630-0017	
Medium – 3 oz			
Otter		630-0020	
Seal	635-0039	630-0021	
Dolphin		630-0022	
Turtle		630-0023	
Medium Heavy – 3.5 oz			
Chipmunk	635-0042	630-0030	636-0113
Rabbit	635-0043	630-0031	636-0114
Fox	635-0044	630-0032	636-0115
Penguin	635-0045	630-0033	636-0116
Monkey	635-0046	630-0034	636-0117
Donkey		630-0035	
Llama		630-0036	
Giraffe		630-0037	



Zoo Pack Elastics

Animal Zoo Packs	Pouches of 15 Zoo Packs	Master Pouches of 50 Zoo Packs	Color Packs of 40 Zoo Packs
Heavy – 4.5 oz			
Kangaroo	635-0050	630-0040	
Bear	635-0051	630-0041	636-0122
Zebra	635-0052	630-0042	
Camel	635-0053	630-0043	
Extra Heavy Weight 6 oz			
Impala	635-0054	630-0050	636-0125
Ram	635-0055	630-0051	636-0126
Moose		630-0052	636-0127
Buffalo		630-0053	636-0128
Light – 8 oz Headgear			
Cougar	635-0058	634-0000	
Leopard	635-0059	634-0001	
Panther	635-0060	634-0002	
Tiger	635-0061	634-0003	
Lion	635-0062	634-0004	
Heavy – 14 oz Headgear			
Walrus	635-0065	634-0012	
Elephant	635-0066		
Whale	635-0067	634-0014	

Zoo Pack Dispenser

- Convenient, clear plexi rack
- Holds 6 Zoo Pack boxes
- Easily mounts on wall or stands on counter top

Dispenser (Zoo Packs not included)

630-0000

